

Microcontrol Laser Therapy Unit.

Salient Features

- Micro controller base for optimum accuracy.
- 80 Clinical protocol Programs.
- 19 User programmable protocol.
- Choice of Two Different Wave length probe.
 - A) IR probe : 830nm/40mW.
 - B) Red Visible probe : 650nm/10mW.
- Back Light LCD Display.
- Continuous and Pulse mode.
- Sensor for LASER output Check.
- Auto time Calculations.
- ABS Shock Proof cabinet.
- Special Compact design for Home visit, clinics, & Hospital.

