

ORDERING INFORMATION

Intelect Mobile Stimulation

2777 - Intelect Mobile Stim

Standard Accessories

10648 - Nylatex Wrap

79967 - 6 x 8cm carbon electrodes (4x)

79970 - 6 x 8cm sponges (4x)

27933 - User Manual on CD

TECHNICAL SPECIFICATIONS

Mains Power: 120-240VAC, 50-60Hz

Weight: 2.3 kg

Dimensions: (L x W x H) 33 x 29 x 16.3cm

Electrical Safety Class: Class 1, Type BF

Safety Tests: IEC 60601-1, IEC 60601-1-2



Intelect Mobile Stim

- 2 channels of stimulation output
- Independent intensity and parameter controls for each channel
- Multiple waveforms on board: (2 Pole IFC, 4 Pole IFC, Galvanic, High Voltage Pulsed Current (HVPC), Microcurrent, Russian, Träbert, Monophasic Triangular, Monophasic Rectangular, Diadynamic, VMS and TENS)
- 15 user-defined memory positions
- Lightweight design, along with battery-powered option and custom designed mobile carrying bag, make transport easy
- Innovative design allows for cart, table top, wall or mobile use





ORDERING INFORMATION

Optional Accessories

- 27478 - Battery Pack
- 27467 - Intelect Mobile Carrying Bag
- 2775 - Therapy System Cart
- 2774 - Vacuum Module w/Therapy System cart
- 2883 - Cart Adapter

WAVEFORMS

| | |
|--|---|
| | Interferential • Four Pole |
| | Interferential • Two Pole (Premodulated) |
| | TENS • Asymmetrical Biphasic • Symmetrical Biphasic |
| | VMS • Continuous • Cycle Timed • Burst |
| | High Voltage Pulsed Current |
| | Microcurrent |
| | Galvanic/Direct Current |
| | Russian |
| | Träbert/Ultra Reiz |
| | Monophasic Triangular |
| | Monophasic Rectangular |
| | Diadynamic |

