

11 COMBINATION THERAPY

INTELECT® MOBILE COMBINATION

Ordering Information

Intellect Mobile Combination

2778 – Intellect Mobile Combo

Standard Accessories

- 27335 - 5 cm² Applicator
- 4248 - Ultrasound Gel
- 10648 - Nylatex® Wrap
- 79967 - 6 x 8cm carbon electrodes (4x)
- 79970 - 6 x 8cm sponges (4x)
- 27979 - User Manual on CD













Optional Accessories

- 27333 - 1 cm² Applicator
- 27334 - 2 cm² Applicator
- 27336 - 10cm² Applicator
- 27478 - Battery Pack
- 27467 - Intellect Mobile Carrying Bag

Technical Specifications

Mains Power: 120-240VAC, 50-60Hz
 Weight: 2.5 kg
 Dims: (L x W x H) 33 x 29 x 25.7 cm
 Electrical Safety Class: Class 1,
 Electrotherapy Type BF, Ultrasound TypeB
 Safety Tests: IEC 60601-1, IEC 60601-1-2,
 IEC 60601-2-5, IEC 60601-2-10

Waveforms

Interferential • Four Pole	
Interferential • Two Pole (Premodulated)	
TENS • Asymmetrical Biphasic • Symmetrical Biphasic	
VMS • Continuous • Cycle Timed • Burst	
High Voltage Pulsed Current	
Microcurrent	
Galvanic/Direct Current	
Russian	
Trabert	
Monophasic Triangular	
Monophasic Rectangular	
Diadynamic	



- ◆ 1 and 3 MHz ultrasound frequencies
- ◆ Pulsed and continuous ultrasound therapy operation (10%, 20%, 50% & 100%)
- ◆ Variable duty cycle outputs of 16Hz, 48Hz, or 100Hz for ultrasound
- ◆ Multiple waveforms on board: (2 Pole IFC, 4 Pole IFC, Galvanic, High Voltage Pulsed Current (HVPC), Microcurrent, Russian, Trabert, Monophasic Triangular, Monophasic Rectangular, Diadynamic, VMS and TENS)
- ◆ 15 user-defined memory positions
- ◆ Lightweight design, along with battery- powered option and custom designed mobile carrying bag, make transport easy
- ◆ Unit uniquely designed for tabletop, wall mount or mobile use
- ◆ Ergonomically designed sound heads for ultrasound therapy in variety of sizes (1 cm², 2 cm², 5 cm², and 10 cm²)
- ◆ All sound heads are interchangeable due to Electronic Signature™ residing on all applicators
- ◆ Head warming feature and head contact control with visual and auditory feedback on all sound heads

