

BIOWAVE HM-801

Microwave Diathermy

Microwave penetrates deeper into tissue, muscle and thereby enhancing local blood circulation in muscle., Subcutaneous organization and others can be warmed evenly without heating skin excessively.

Microwave therapy can be applied to treatment portion of a patient with in dress just by turning radiator to treatment portion. Biowave HM-801 has versatile concentrated heat radiation system and the success of therapy including thermal effect is indisputable.

Features

- Easy controlling power on/off and strength and simple installing circular or oblong irradiator
- Any angle and direction for any treatment part
- Operation starts only in zero intensity
- Overheating(80°C) is prevented with magnetron protective circuit.

Specification

- Power : AC 220V, 60Hz
- Output Power : 0 ~ 240 W
- Output Frequency : 2450 MHz
- Wavelength : 124.4 nm
- Output Mode
 - 1) Continuous
 - 2) Intermittent
- Weight : 22 kg
- Dimension : 430(W) x 315(D) x 760(H) mm

