

TTET-200 and TTFT-200 Traction Tables



The TTET-200 is a two-section, 71 cm wide table with an adjustable height from 48 cm - 96 cm (19" - 38") and a working load of 181 kg (400 lbs). Optional retractable casters are available. Traction unit sold separately. Available in all upholstery colors shown in chart.



The TTFT-200 is a two-section, 71 cm wide table with a fixed height 81 cm (32") from top to floor and a working load of up to 159 kg (350 lbs). Traction unit sold separately. Available in all upholstery colors shown in chart.



ORDERING INFORMATION

Economy Hi-Lo Traction Table

6870 – TTET-200

Optional Accessories

66898 – Foot Switch

86016 – Retractable Casters

TECHNICAL SPECIFICATIONS

Mains Power: 230 V, 50/60 Hz;
120 V, 60 Hz

Current: 230 V - 1.0 A; 120 V - 1.6 A

Power Input: 230 V - 230 W;
120 V - 200 W

Weight: 64 kg (140 lb)

Dimensions: (L x W x H) 196 x 71 x
48 to 96 cm (77" x 28" x 20.5" to 41")

Electrical Safety Class: Class 1, Type B

Economy Fixed Height Traction Table






6840 – TTFT-200

UPHOLSTERY*

Standard Colors

	Black (002)
	Imperial Blue (001)
	Mocha (003)
	Burgundy (004)
	Forest (005)
	Slate (006)
	Midnight Gray (033)
	Dove (007)
	Stone Blue (061)

Designer Colors

	Seafoam (026)
	Sandstone (027)
	Terra Cotta (029)
	Blue Ridge (032)
	Rosewood (061)

Rhinohyde Colors

	Northwoods (052)
	Cobalt (053)
	Port (056)
	Plum (059)

*See page 6 for more information.